ACUMED®

ACUMED CLAVICLE SCREW SYSTEM

MediMark® Europe Sarl.
11 rue Emile ZOLA. BP 2332
38033 GRENOBLE CEDEX 2
FRANCE
+33 4 76 86 43 22

Acumed®
5885 NW Cornelius Pass Road
Hillsboro, OR 97124-9432
+1 (503) 627-9957
www.acumed.net

PKGI-53-B
Effective 11/2009

0473

* P K G I - 5 3 - B *

ACUMED SKRUESYSTEM TIL CLAVICULA

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ACUMED CLAVICLE SCREW SYSTEM

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

DESCRIPTION: The Acumed Clavicle Screw System is designed to provide fixation to repair an acute fracture, mal-union or non-union of the clavicle, during bone healing.

INFORMATION FOR USE: The surgeon must select the type and size implant that best meets the patient's surgical needs.

INDICATIONS: The Acumed Clavicle Screw System is indicated for fusions, fractures, or osteotomies of the clavicle.

CONTRAINDICATIONS: Contraindications for this system are active or latent infection, sepsis, osteoporosis, insufficient quantity or quality of bone/soft tissue, and material sensitivity. If sensitivity is suspected, tests are performed prior to implantation. Patients who are unwilling or incapable of following postoperative care instructions are contraindicated for this device. This device is not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical or lumbar spine.

WARNINGS: For safe effective use of this implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for this device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, non-union, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result

of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, non-union, or incomplete healing, and the possibility of nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail.

PRECAUTIONS: An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Instruments should be inspected for wear or damage prior to usage. Protect implants against scratching and nicking, as such stress concentrations can lead to failure. Particular care should be paid to hex drivers, drill bits and instruments used for implant insertion.

ADVERSE EFFECTS: Possible adverse effects are pain, discomfort, or abnormal sensations due to the presence of an implant. Implant fracture, migration and/or loosening may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Metal sensitivity or histological or allergic reaction resulting from implantation of a foreign material may occur. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing may result from the presence of an implant or due to surgical trauma.

STERILITY: The product is provided sterile and non-sterile. Sterile device was exposed to a minimum dose of 25.0-kGy gamma irradiation. Resterilisation may only be performed if the original sterile package has been opened in error using one of the following methods.

Sterilization Methods

Gravity Displacement Autoclave	270° F (132° C) for 40 minutes
Gravity Displacement Autoclave Dry Time	70 minutes
Pre-vacuum Autoclave	270° F (132° C) for 14 minutes
Pre-vacuum Autoclave Dry Time	70 minutes

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

CAUTION: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital.

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CONTRAINDICATIONS: Contraindications for this system are active or latent infection, sepsis, osteoporosis, insufficient quantity or quality of bone/soft tissue, and material sensitivity. If sensitivity is suspected, tests are performed prior to implantation. Patients who are unwilling or incapable of following postoperative care instructions are contraindicated for this device. This device is not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, or lumbar spine.

WARNINGS: For safe effective use of this implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for this device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result

of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing, and the possibility of nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail.

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STERILITY: The product is provided sterile and non-sterile. Sterile device was exposed to a minimum dose of 25.0-kGy gamma irradiation. Resterilization may only be performed if the original sterile package has been opened in error using one of the following methods.

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